

Mobile N2O Support System

Product Code: **N20MSS**



Specifications:

NITROUS OXIDE

Twin "E" Size Cylinder Trolley

5 Star Trolley Base for Increased Stability

Manual Twin Cylinder Change Over System

British Standard Gas Indexed Coloured Self Sealing Valve

Complete With:

2x R24105N N2O Pin-Index Compact Regulators

Weights/Dimensions:

Height: 940mm

Width: 600mm

Depth: 600mm

Weight: 10kg (approx)

Oxylitre Holdings Limited offer the most comprehensive range of Pipeline and Hose Accessories. Our Pipeline and Hose Accessories cover all Medical Gases suitable. Our connections are suitable for use with all flexible medical hose. Our extensive range of gas specific pipeline connections and accessories cover both British and European Standard connections.

Oxylitre's range of Pipeline and Hose Accessories are dependable, extremely easy to maintain with a long service life and a **7 year** warranty.