

Compact Regulator Pin Index and Self Sealing Valve Oxygen

Product Codes: **R14015-S & R24015-S**



Specifications:

OXYGEN

Stainless Steel Body

Pin Index Cylinder Fitting

138 bar maximum inlet pressure (R14015-S)

200 bar maximum inlet pressure (R24015-S)

Pre-Set 4 bar nominal output pressure

High Pressure Luminous Cylinder Content Gauge

BS Gas Indexed Coloured Self Sealing Valve

Complete with High Pressure Sensor

7 Year Manufacturers Warranty

High Pressure Sensor Specification:

Range: 0-200 Bar

Power Supply: 8-30 VDC

Output: 1-5 V

Accuracy: 0.5% (+/- 0.5% FS at 25°C)

IP Rating: IP67

4 Pin Connecting Wire, M12x1 Nut

Wire Length: 2.5m (custom lengths available)

Weights/Dimensions:

Weight: 0.83 kg

Height: 150mm

Width: 124mm

Depth: 62mm

Accessories:

SK40 - 5 Year Service Pack

Oxylitre' range of Compact Medical Regulators, have been designed with the end user in mind. Our 138 bar rated Compact Range of regulators are simple to use, reliable and extremely functional, Oxylitre' 138 bar Compact Range Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Compact Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.